

# TRANSDUCTION OF ANTIGEN-SPECIFIC CD8<sup>+</sup> T CELLS WITH RETROVIRAL VECTORS

Francesco Marangoni January 13, 2010

An easy and flexible protocol to transduce CD8<sup>+</sup> T cells and differentiate them into central or effector memory cells.

---

## MATERIALS

ACK buffer to lyse red blood cells (NH<sub>4</sub>Cl 150mM, KHCO<sub>3</sub> 10mM Na<sub>2</sub>EDTA 0.1mM)

T cell medium: RPMI 10% FCS HEPES (1:66) Sodium Pyruvate (1:100) GlutaMax (1:100) β-mercaptoethanol (1:1000)

Cognate peptide (e.g. OT-I or CL4 peptide)

Polybrene solution 4 mg/ml

If cells are from Balb/c, you need IL-12 (R&D 419-ML)

To generate effector memory cells, you need IL-2 (R&D 1150-ML)

To generate central memory cells, you need IL-15 (R&D 447-ML)

---

## METHOD

### Day 1:

#### 1. Culture and stimulate CD8<sup>+</sup> T cells:

- a. Sacrifice mice, take spleen and LN (axillary, brachial, inguinal).
- b. Smash organs in sterile conditions. Put them in a petri dish and keep wet with PBS. LNs can be ground using 2 frosted and ethanol-wiped glass slides (be careful not to break them). Filter the suspension (40µm filter) and thoroughly wash the Petri dish. Centrifuge the cells.
- c. Lyse red blood cells by resuspension in 1ml ACK for 5' at room temperature. Block lysis by adding 9ml PBS.
- d. Resuspend cells in 10ml T cell medium + 10 µM of cognate peptide.
- e. **PRO TIP:** If cells come from B6 mice, APC produce enough IL-12. If they come from Balb/c, this is not the case. Supplement IL-12 20ng/ml.
- f. Put in a 10-cm Petri dish and culture overnight.

### Day 2:

#### 1. First transduction:

- a. Thaw retroviral supernatant (10ml per sacrificed mouse) and add polybrene 4 µg/ml. If cells are of Balb/c origin, also add 20ng/ml IL12.
- b. Harvest the T cells stimulated yesterday and spin them down at 500g for 5'.
- c. Resuspend the cells in the viral supernatant prepared in step a, and aliquot them 1ml/well of a 24-well plate.
- d. Spinfect 1000g 90' at 32C.
- e. After spinfection, add to each well 1ml of T cell medium (+ 20ng/ml IL-12 if cells are from Balb/c mice).
- f. Incubate overnight.

### Day 3:

#### 1. Second transduction:

- a. Repeat Day 2 steps 1a to 1d.

- b. After spinfection, accurately harvest cells, which are stuck to the bottom of the wells, and transfer them in a 15ml Falcon tube. Spin 500g 5'.
- c. Resuspend cells in T cell medium with 20ng/ml IL-2, to generate effector memory cells.
- d. Plate cells in Petri dishes at the concentration of  $1 \times 10^6$ /ml.  
**PRO TIP:** Substituting IL-2 with 20ng/ml IL-15 will generate central memory cells. In this case maintaining cell concentration to at least  $1 \times 10^6$ /ml is critical, as IL-15 is trans-presented by APC still present in the culture.

*Day 5:*

1. Select transduced cells. This can occur in three ways, depending on the selection gene encoded by the retroviral vectors:
  - i. Retro encodes puromycin resistance: apply puromycin 8  $\mu$ g/ml;
  - ii. Retro encodes  $\Delta$ LNGFR: do immunomagnetic selection for  $\Delta$ LNGFR<sup>+</sup> cells, or
  - iii. Retro encodes fluorescent proteins: FACS sort cells.
2. Culture cells in T cell medium with 20ng/ml IL-2 or 20ng/ml IL15.

*Day 6 or 7:*

1. Use cells in experiments.